### Faster clearance, fewer treatments





















PicoSure. The world's first picosecond laser.

about a PicoSure treatment that's right for you.

Learn more about PicoSure at www.picosure.com




PicoSure 755 is US FDA cleared to treat tattoos and pigmented lesions in skin types I-VI. PicoSure 532 is US FDA cleared to treat tattoos in skin types I-III. Results may vary.

Cover model used for illustrative purposes. Results may vary.

©2016 Cynosure, Inc. All rights reserved. Cynosure and PicoSure are registered trademarks and PressureWave is a trademark of Cynosure, Inc.









### Ready to lose that tattoo? Do it better and faster with PicoSure.®

PicoSure is the world's first picosecond laser and the only advanced tattoo removal laser with PressureWave $^{\text{TM}}$  technology, which works with your body to remove unwanted ink.

#### Selective "shattering" with PressureWave™

PicoSure delivers laser energy that specifically targets tattoo ink, resulting in the successful shattering of just the ink particles without harming the surrounding tissue. PressureWave is the secret behind tattoo treatments with PicoSure.

PicoSure works effectively on:


- Multicolored tattoos, including blues and greens
- · Hard-to-remove black ink
- Stubborn, previously treated tattoos


Best of all, PicoSure tattoo removal often requires fewer treatments than traditional laser tattoo removal.

# How does a PicoSure treatment work?

During treatment with PicoSure, ultra-short (trillionths of a second) pulses target ink particles with PressureWave technology—a pressure so powerful, the ink shatters into tiny dust-like particles. These microscopic particles are more rapidly absorbed and eliminated naturally by the body.

#### PressureWave technology in action







Electron Microscopy Courtesy of H.B. Jalian, MD

Unlike traditional tattoo lasers, which rely on heat to destroy tattoo ink, PicoSure delivers a photomechanical impact called PressureWave that effectively shatters ink into tiny particles.

This can mean better clearance of the tattoo and fewer treatments compared to conventional laser tattoo removal.

# Lose the ink, faster, with PicoSure

- Fast treatment times of typically 10 minutes
- Number of treatments needed is dependent on variables such as ink characteristics, size of tattoo, location on body, and lifestyle habits\*
- Treatments typically occur every 6-8 weeks\*
- \* Your treatment provider will develop a treatment plan to meet your individual needs.

Patients and treatment providers have reported better clearance (disappearance) of tattoos with improved recovery time. So you can expect excellent results with PicoSure.

PressureWave your tattoo away with PicoSure

Let PicoSure work with your body to clear unwanted ink colors, with breakthrough tattoo removal.

PicoSure works with your body to remove tattoos—better and in fewer treatments.

